



# 2021

## Wisconsin Milk Proficiency Testing Drug Residue Set

### *Instructions for Handling Test Materials*

#### Result Forms

- Download the fillable PDF [Drug Residue Result Form](#). Adobe Acrobat or Adobe Acrobat Reader will be needed to complete this form. The latest version of Adobe Acrobat Reader can be downloaded for free at [Adobe](#).
- Full details regarding entering and reporting results are described in each result form and the [Proficiency Testing Forms Instructions](#).
- Submit result forms by **April 9, 2021**.

#### Receipt and Storage

- Verify samples 15-23 are present.

Test Set Contents		
Sample(s)	Vial Color	Test
15 – 22	Green	Drug Residue (Inhibitory Substances)
23	Clear	Temperature Control

- Record temperature of TC (#23) on test result form (page 1).
- If any samples are frozen, record sample number on test result form (page 1). Frozen samples may be thawed and tested.
- If any samples are leaking, please contact [DATCPBLSProficiencyTesting@wisconsin.gov](mailto:DATCPBLSProficiencyTesting@wisconsin.gov) or (608) 224-4831.
- Keep all samples refrigerated (0 – 4.5 °C).
- Analyze all samples within one day of receipt.

#### Analysis

- Test materials should be analyzed in accordance with the most current FDA 2400 series.
- Keep all samples within 0 – 4.5 °C during analysis.
- Each analyst should perform test controls independently.
- Samples should not be analyzed in duplicate for any test.

#### Shipping Supplies

- Return shipper to:  
  
WDATCP Bureau of Laboratory Services  
2601 Agriculture Drive, Suite 150  
Madison, WI 53718
- The follow items must be returned with the shipper: icepack, styrofoam square, polycell floater vial rack, and pegboard.
- **Please note that if shippers are not returned by April 19, 2021 a \$50 replacement fee will be added to the billing invoice.**