



Wisconsin Milk Proficiency Testing Drug Residue Set

Instructions for Handling Test Materials

Result Forms

- Download the fillable PDF **Drug Residue Result Form** from the [Milk Proficiency Testing Program](#) website. Adobe Acrobat or Adobe Acrobat Reader will be needed to complete this form. The latest version of Adobe Acrobat Reader can be downloaded for free at [Adobe](#).
- Full details regarding entering and reporting results are described in each result form and the **Proficiency Testing Forms Instructions** (located on the [Milk Proficiency Testing Program](#) website).
- See website for results submission deadline.

Receipt and Storage

- Verify samples 15 – 23 are present.

Test Set Contents		
Sample(s)	Vial Color	Test
15 – 22	Clear	Drug Residue (Inhibitory Substances)
23	Colored	Temperature Control

- Record temperature of TC (#23) prior to testing on test result form (page 1).
- If samples are frozen, record sample number on test result form (page 1). Frozen samples may be thawed and tested.
- If any samples are leaking, please contact DATCPBLSProficiencyTesting@wisconsin.gov or (608) 224-4831.
- Keep all samples refrigerated (0 – 4.5 °C).
- Analyze all samples within one day of receipt.**

Analysis

- Test materials should be analyzed in accordance with the most current NCIMS 2400 series.
- Keep all samples within 0 – 4.5 °C during analysis.
- Each analyst should perform test controls independently.
- Samples should not be analyzed in duplicate for any test.

Shipping Supplies

- Return shipper to:

WDATCP Bureau of Laboratory Services
2601 Agriculture Drive, Suite 150
Madison, WI 53718
- The follow items must be returned with the shipper: icepacks, polycell floater vial rack, and pegboard(s).
- Please note that if shippers are not returned by May 1st a \$50 replacement fee will be added to your billing invoice.**